

Instructions for Use

Cleaning & Disinfection: Koda XT

Compatible Cleaners & Disinfectants*

Liquids

- Diversey™ Oxivir® Tb Hospital Grade Disinfectant
- Diversey™ Avert® Disinfectant Cleaner
- Ecolab® Oxycide Daily Disinfectant Cleaner
- ChemiSwiss SUclean™
- ClearSpace powered by PreventX™
- Hydrogen Peroxide (3% Solution)
- Isopropyl Alcohol

Avoid using a dry cloth, scrapers, or squeegees on 3form Koda XT. For more information on standard recommended cleaning procedures, refer to the Koda XT Spec Sheet available on www.3-form.com

Koda XT is an exterior grade resin designed to withstand the demands of both exterior environmental conditions and the rigorous cleaning and disinfection requirements in Healthcare.

Koda XT provides excellent chemical resistance which allows for effective cleaning and disinfection of environmental surfaces to reduce microbe transmission. Refer to the compatible cleaners and disinfectants to ensure the proper disinfection and minimal wear.

